


Respireo
NINA



Silent and intuitive

UNIVEO Δ headgear

The first innovative headgear designed to help the patient

Quick Adjust, the Univeo Delta headgear scaling feature (patented)



Landmarks on the lower straps

Quick Adjust feature with landmarks

Easy to handle and set

Soft & seamless fabric

Two-tone

Facilitates the installation and adjustment of the headgear

► **helps to adjust the mask and promotes adherence**

A generation of
innovative masks

Respireo
NINA

Univeo Delta headgear

With **Quick Adjust** scaling function
A headgear designed to help patients
Available in S and M/L



Frame

Minimalist and light
Optimal coupling
ergonomics

Quiet Flow

The function that ensures
optimal silence
Only 12.3 ± 3.2 dB(A)
sound pressure level

Secure Clip

The mistake-proof system

Cushion

Optimal sealing and comfort
Available in S, M and L sizes

Ergonomic hooks

Quick assembly and
disassembly of the headgear



Made in Italy

Respireo NINA

2 headgears
3 cushions
1 frame



Product Reference	Description	Size	Code for order
472367	Respireo NINA	S	KM229116
472368	Respireo NINA	M	KM229117
472369	Respireo NINA	L	KM229118

Spare parts:

472478	Respireo NINA cushion	S	KM229122
472479	Respireo NINA cushion	M	KM229123
472480	Respireo NINA cushion	L	KM229124
472481	Univeo Delta headgear	S	KM229125
472482	Univeo Delta headgear	M/L	KM229126
472483	Swivel connector	-	KM229130



Watch
the video

Manufactured by

AIR LIQUIDE MEDICAL SYSTEMS S.r.l.
Via dei Prati, 62
25073 - Bovezzo (BS) ITALY

Distributor contact

AIR LIQUIDE MEDICAL SYSTEMS S.A.
Parc de Haute Technologie
6, rue Georges Besse
92182 Antony Cedex
France

www.device.airliquidehealthcare.com
Mail: dl-frtaema-src-export@airliquide.com

Hotline ALMS: +33 179 51 70 01

Air Liquide
HEALTHCARE

Respireo NINA is a reusable nasal mask with calibrated exhalation orifices (Vented) intended to be used at home or in a hospital setting by a single adult patient weighing over 30 Kg, who have been prescribed with a non-invasive positive pressure ventilation therapy (NIPPV), e.g. continuous positive airway pressure therapy (CPAP) or Bilevel positive airway pressure therapy (BiPAP).

Medical device class IIa CE0051. Please read the user manual.